

**BoneScalpel™**  
**Sterilization Instructions\***



**With Sterilization Tray**

	132°C (270°F)	134-137°C (274-279°F)
Configuration	Items placed in Misonix Sterilization Tray E-SYSTRAY <sup>1</sup> or E-SYSTRAY-2 <sup>1</sup>	Items placed in Misonix Sterilization Tray E-SYSTRAY <sup>2</sup> or E-SYSTRAY-2 <sup>2</sup>
Cycle	Prevacuum	Prevacuum
Preconditioning Pulses	4	4
Minimum Exposure Time	4 minutes*	3 minutes*
Minimum Dry Time	30 minutes	30 minutes

**Wrapped, No Tray**

	132°C (270°F)	134-137°C (274-279°F)	134°C/273°F
Configuration	Items wrapped, NO TRAY <sup>3</sup>	Items wrapped, NO TRAY <sup>3</sup>	Items wrapped, NO TRAY <sup>3</sup>
Cycle	Prevacuum	Prevacuum	Gravity
Preconditioning Pulses	4	4	None
Minimum Exposure	4 minutes*	3 minutes*	20 minutes
Minimum Dry Time	45 minutes	30 minutes	5 minutes

<sup>1</sup>Tray wrapped in Kimberly Clark KC300 KIMGUARD sterilization wrap.

<sup>2</sup>Tray wrapped in Kimberly Clark KC300 or KC600 KIMGUARD sterilization wrap.

<sup>3</sup>Items wrapped in Kimberly Clark KC300 or KC600 KIMGUARD sterilization wrap.

\*Exposure time can be increased up to a maximum of 18 minutes to comply with local requirements and/or recommendations of the World Health Organization (WHO), Robert Koch Institute (RKI), etc. Misonix Inc reusable medical devices are able to sustain such sterilization cycles.

WI-10-210 Rev A