

SonicOne® O.R Ultrasonic Surgical Debridement System

*Sterilization Instructions**

With Sterilization Tray

	132°C (270°F)	134-137°C (274-279°F)
Configuration	Items placed in Misonix Sterilization Tray E-SYSTRAY ¹ or E-SYSTRAY-2 ¹	Items placed in Misonix Sterilization Tray E-SYSTRAY ² or E-SYSTRAY-2 ²
Cycle	Prevacuum	Prevacuum
Preconditioning Pulses	4	4
Minimum Exposure Time	4 minutes*	3 minutes*
Minimum Dry Time	30 minutes	30 minutes

Wrapped, No Tray

	132°C (270°F)	134-137°C (274-279°F)	134°C/273°F
Configuration	Items wrapped, NO TRAY ³	Items wrapped, NO TRAY ³	Items wrapped, NO TRAY ³
Cycle	Prevacuum	Prevacuum	Gravity
Preconditioning Pulses	4	4	None
Minimum Exposure	4 minutes*	3 minutes*	20 minutes
Minimum Dry Time	45 minutes	30 minutes	5 minutes

¹Tray wrapped in Kimberly Clark KC300 KIMGUARD sterilization wrap.

²Tray wrapped in Kimberly Clark KC300 or KC600 KIMGUARD sterilization wrap.

³Items wrapped in Kimberly Clark KC300 or KC600 KIMGUARD sterilization wrap.

*Exposure time can be increased up to a maximum of 18 minutes to comply with local requirements and/or recommendations of the World Health Organization (WHO), Robert Koch Institute (RKI), etc. Misonix Inc reusable medical devices are able to sustain such sterilization cycles.